

# RASGRP3 Rabbit pAb

货号: **AYP25231**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>WB:</b> Homo sapiens
应用	<a href="#">WB</a> <a href="#">IF/ICC</a>
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IF/ICC:</b> 1:50 - 1:200
理论分子量	78kDa
实测分子量	78kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BT-474,A-549,Mouse heart,Mouse lung,Mouse thymus
细胞定位	perinuclear region of cytoplasm,plasma membrane
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 540-689 of human RAS GRP3 (NP_056191.1).
序列	VLACRRFARAPSLSSGHGSLPGSPSLPPAQDEVFEFPGVTAGHRDLDSRAITLVTGSSRKISVRLQRATTSQATQTEPVWS EAGWGDSGSHTFPKMKSKFHDKAAKDKGFAKWENKPRVHAGVDVDRGTEFELDQDEGEETRQDGEDG

## 靶点信息

研究背景	The protein encoded by this gene is a guanine nucleotide exchange factor that activates the oncogenes H RAS and RAP1A. Defects in this gene have been associated with systemic lupus erythematosus and several cancers.
基因ID	25780

基因名	RASGRP3
Swiss	Q8IV61
别名	RASGRP3;GRP3

## 产品验证



Western blot analysis of RASGRP3 expressed in Mouse heart,BT-474,A-549 using RASGRP3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)