

# STXBP3 Rabbit pAb

货号: **AYP25305**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	67kDa
实测分子量	68kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,U-251MG,Jurkat,Mouse kidney,Mouse lung,Mouse spleen,Mouse brain,Rat brain
细胞定位	Cell membrane,Cytoplasm,cytosol
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 343-592 of human STXB P3 (NP_009200.2).
序列	KQVVHLNLAEDCMNKFKNLIEKLCCKTEQDLALGTD AEGQKVKDSMRVLLPVLLNKNHDNCDKIRAILLYIFSINGTTEENLD RLIQNVKIENESDMIRNWSYLGVPVIPSQQGKPLRKDRSAEETFQLSRWTPFIKDIMEDAIDNRLDSKEWPYCSQCPAVW NGSGAVSARQKPRANYLEDKNGSKLIVFVIGGITYSEVRCAYEVSQAHKSCVEIIGSTHVLTPKLLDDIKMLNPKDKVSLIKDE

## 靶点信息

研究背景	Together with STX4 and VAMP2, may play a role in insulin-dependent movement of GLUT4 and in docking /fusion of intracellular GLUT4-containing vesicles with the cell surface in adipocytes.
基因ID	6814

基因名	STXBP3
Swiss	O00186
别名	STXBP3;MUNC18-3;MUNC18C;PSP;UNC-18C

## 产品验证



Western blot analysis of STXBP3 expressed in Mouse kidney,Rat brain,U-87MG using STXBP3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)