

— ABLYBIO, Help Your Research



# PHF21B Rabbit pAb

货号: **AYP25322**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	51kDa/53kDa/57kDa
实测分子量	50kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Mouse spinal cord,Rat brain
细胞定位	
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 322-531 of human PHF21B (NP_612424.1).
------	--

## 靶点信息

研究背景	Predicted to enable metal ion binding activity.
------	---

基因ID	112885
基因名	PHF21B
Swiss	Q96EK2 ( <a href="https://www.uniprot.org/uniprotkb/Q96EK2/entry">https://www.uniprot.org/uniprotkb/Q96EK2/entry</a> )
别名	PHF21B,BHC80L,PHF4,PHF21B Rabbit pAb,KIAA1661

## 产品验证



Western blot analysis of PHF21B expressed in Mouse brain, Mouse spinal cord, Rat brain using PHF21B Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)