

PHF21B Rabbit pAb

货号: AYP25322

产品信息

| | |
|-------|--|
| 反应 | Human,Mouse,Rat |
| 宿主 | Rabbit |
| 克隆性 | Polyclonal |
| 预测反应 | |
| 应用 | WB |
| 推荐浓度 | WB: 1:500 - 1:2000 |
| 理论分子量 | 51kDa/53kDa/57kDa |
| 实测分子量 | 50kDa |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | Mouse brain,Mouse spinal cord,Rat brain |
| 细胞定位 | |
| 纯化 | Affinity purification |

抗原信息


| | |
|------|--|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 322-531 of human PHF21B (NP_612424.1). |
| 序列 | KRLASNYLNNPLFLTARANEDPCWKNEITHDEHCAACKRGANLQPCGTCPGAYHLSCLEPPLKTAPKGVWVCPRCQQKAL KKDEGVPWTGMLAIVHSYVTHKTVKEEEKQKLLQRGSELQNEHQLEERDRRLASAVQKCLELKTSLARQRGTQSSSLDR LRALLRLIQGEQLLQVTMTTTPAPLLAGPWTKPSVAATHPTVQHPQGHN |

靶点信息

| | |
|------|---|
| 研究背景 | Predicted to enable metal ion binding activity. |
| 基因ID | 112885 |
| 基因名 | PHF21B |

| | |
|--------------|--------------------|
| Swiss | Q96EK2 |
| 别名 | PHF21B;BHC80L;PHF4 |

产品验证



Western blot analysis of PHF21B expressed in Mouse brain, Mouse spinal cord, Rat brain using PHF21B Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn