

HLF Rabbit pAb

货号: **AYP25354**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	23kDa/33kDa
实测分子量	36kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver,Mouse brain
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-295 of human HLF (NP_002117.1).
序列	MEKMSRPLPLNPTFIPPPYGVLRSLLLENPLKLPLHHEDAFSKDKDKEKKLDDESNSPTVPQSAFLGPTLWDKTLPYDGDTFQLEYMDLEEFLESENGIPSPSQHDHSPHPPGLQPASSAAPSVM DLSSRASAPLHPGIPSPNCMQSPIRPGQLLPANRNTPS PIDPDTIQVPVGYEPDPADLALSSIPGQEMFDPKRKRKFSEELKPQPMIKKARKVFIPDDLKDDKYWARRRRKNNMAAKRSR DARRLKENQIAIRASFLEKENSALRQEVADLRKELGKCKNILAKYEARHGPL

靶点信息

研究背景	This gene encodes a member of the proline and acidic-rich (PAR) protein family, a subset of the bZIP transcription factors. The encoded protein forms homodimers or heterodimers with other PAR family members and binds sequence-specific promoter elements to activate transcription. Chromosomal translocations fusing portions of this gene with the E2A gene cause a subset of childhood B-lineage acute lymphoid leukemias. Alternatively spliced transcript variants have been described, but their biological validity has not been determined.
------	---

基因ID	3131
基因名	HLF
Swiss	Q16534
别名	HLF

产品验证



Western blot analysis of HLF expressed in Mouse liver, Mouse brain using HLF Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn