

EYA3 Rabbit pAb

货号: **AYP25370**

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	48kDa/56kDa/57kDa/58kDa/62kDa
实测分子量	57kDa/58kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,Rat eye
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human EYA3 (NP_001981.2).
序列	MEEEQDLPEQPVKKAKMQESGEQTISQVSNPDVSDQKPETSSLASNLPMSEEIMTCTDYIPRSSNDYTSQMYSAPYAHIL SVPVSETAYPGQTQYQTLQQTQPYAVYPQATQTYGLPPFGALWPGMKPESGLIQTSPSPQHSVLTCTTGLTTSQPSPAHY SYPIQASSTNASLISTSSTIANIPAAAVASISNQDYPTYTILGQNQYQACYPSSSFGVTGQTNSDAESTTLAATTYQSEKPSV MAPAPAAQRLSSGDPSTSPSLSQTTPSKDTDDQSRKNMNTSKNRGKRKADATSSQDS

靶点信息

研究背景	This gene encodes a member of the eyes absent (EYA) family of proteins. The encoded protein may act as a transcriptional activator and have a role during development. It can act as a mediator of chemoresistance and cell survival in Ewing sarcoma cells, where this gene is up-regulated via a micro-RNA that binds to the 3' UTR of the transcript. A similar protein in mice acts as a transcriptional activator. Alternative splicing of this gene results in multiple transcript variants.
------	--

基因ID	2140
基因名	EYA3
Swiss	Q99504
别名	EYA3

产品验证



Western blot analysis of EYA3 expressed in HT-29,Rat eye using EYA3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn