

— ABLYBIO, Help Your Research



VEGFR3/FLT4 Rabbit pAb

货号: **AYP25381**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IF/ICC
推荐浓度	IF/ICC: 1:50 - 1:200
理论分子量	93kDa/145kDa/152kDa
实测分子量	178kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat lung
细胞定位	Cell membrane,Cytoplasm,Nucleus,Secreted,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 25-776 of human VEGFR 3/FLT4 (NP_002011.2).
------	---

靶点信息

研究背景	This gene encodes a tyrosine kinase receptor for vascular endothelial growth factors C and D. The protein is thought to be involved in lymphangiogenesis and maintenance of the lymphatic endothelium. Mutations in this gene cause hereditary lymphedema type IA.
基因ID	2324
基因名	FLT4
Swiss	P35916 (https://www.uniprot.org/uniprotkb/P35916/entry)
别名	FLT4,FLT-4,FLT41,LMPH1A,PCL,VEGFR-3,VEGFR3,VEGFR3/FLT4 Rabbit pAb,Fms-like tyrosine kinase 4,Tyrosine-protein kinase receptor FLT4

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)