

— ABLYBIO, Help Your Research



Aurora A (Phospho-Thr288) Antibody

货号: **AYP4074**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	45kDa
实测分子量	46kD
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,HT-29
细胞定位	Cytoplasm,centriole,centrosome,cilium basal body,cytoskeleton,microtubule organizing center,spindle pole
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human Aurora A (Phospho-Thr288).
------	---

靶点信息

研究背景	The protein encoded by this gene is a cell cycle-regulated kinase that appears to be involved in microtubule formation and/or stabilization at the spindle pole during chromosome segregation. The encoded protein is found at the centrosome in interphase cells and at the spindle poles in mitosis. This gene may play a role in tumor development and progression. A processed pseudogene of this gene has been found on chromosome 1, and an unprocessed pseudogene has been found on chromosome 10. Multiple transcript variants encoding the same protein have been found for this gene.
基因ID	6790
基因名	AURKA
Swiss	O14965 (https://www.uniprot.org/uniprotkb/O14965/entry)
别名	AURKA,AIK,ARK1,AURA,BTAK,PPP1R47,STK15,STK6,STK7,Aurora A (Phospho-Thr288) Antibody,Aurora 2,Aurora/IPL1-related kinase 1,Breast tumor-amplified kinase,lpl1- and aurora-related kinase 1,Serine/threonine-protein kinase 15,Serine/threonine-protein kinase 6,Serine/threonine-protein kinase Ayk1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)