

— ABLYBIO, Help Your Research



Lck (Phospho-Tyr505) Antibody

货号: **AYP4091**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	40kDa/58kDa/61kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,NIH/3T3
细胞定位	Cell membrane,Cytoplasm,Cytoplasmic side,Lipid-anchor
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human Lck (Phospho-Tyr505).
------	--

靶点信息

研究背景	This gene is a member of the Src family of protein tyrosine kinases (PTKs). The encoded protein is a key signaling molecule in the selection and maturation of developing T-cells. It contains N-terminal sites for myristylation and palmitoylation, a PTK domain, and SH2 and SH3 domains which are involved in mediating protein-protein interactions with phosphotyrosine-containing and proline-rich motifs, respectively. The protein localizes to the plasma membrane and pericentrosomal vesicles, and binds to cell surface receptors, including CD4 and CD8, and other signaling molecules. Multiple alternatively spliced variants encoding different isoforms have been described.
基因ID	3932
基因名	LCK
Swiss	P06239 (https://www.uniprot.org/uniprotkb/P06239/entry)
别名	LCK,IMD22,LSK,YT16,p56lck,pp58lck,Lck (Phospho-Tyr505) Antibody,Leukocyte C-terminal Src kinase,Lymphocyte cell-specific protein-tyrosine kinase,Protein YT16,Proto-oncogene Lck,T cell-specific protein-tyrosine kinase,p56-LCK,Lck (Phospho-Tyr505)

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)