

— ABLYBIO, Help Your Research



# NMDAR1 (Phospho-Ser890) Antibody

货号: **AYP4120**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200 <b>IF/ICC:</b> 1:50 - 1:200
理论分子量	99-107kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Rat brain
细胞定位	Cell junction,Cell membrane,Multi-pass membrane protein,postsynaptic cell membrane,postsynaptic density,synapse
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human NMDAR1 (Phospho-Ser890).
------	---

## 靶点信息

研究背景	The protein encoded by this gene is a critical subunit of N-methyl-D-aspartate receptors, members of the glutamate receptor channel superfamily which are heteromeric protein complexes with multiple subunits arranged to form a ligand-gated ion channel. These subunits play a key role in the plasticity of synapses, which is believed to underlie memory and learning. Cell-specific factors are thought to control expression of different isoforms, possibly contributing to the functional diversity of the subunits. Alternatively spliced transcript variants have been described.
基因ID	2902
基因名	GRIN1
Swiss	Q05586 ( <a href="https://www.uniprot.org/uniprotkb/Q05586/entry">https://www.uniprot.org/uniprotkb/Q05586/entry</a> )
别名	GRIN1,GluN1,MRD8,NMD-R1,NMDA1,NMDAR1,NR1,NMDA 1,NMDAR1 (Phospho-Ser890) Antibody,Glutamate [NMDA] receptor subunit zeta-1,N-methyl-D-aspartate receptor subunit NR1,NMDAR1 (Phospho-Ser890)

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)