

— ABLYBIO, Help Your Research



# GATA2 (Phospho-Ser401) Antibody

货号: **AYP4157**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	49kDa/50kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	Nucleus
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human GATA2 (Phospho-Ser401).
------	--

## 靶点信息

研究背景	This gene encodes a member of the GATA family of zinc-finger transcription factors that are named for the consensus nucleotide sequence they bind in the promoter regions of target genes. The encoded protein plays an essential role in regulating transcription of genes involved in the development and proliferation of hematopoietic and endocrine cell lineages. Alternative splicing results in multiple transcript variants.
基因ID	2624
基因名	GATA2
Swiss	P23769 ( <a href="https://www.uniprot.org/uniprotkb/P23769/entry">https://www.uniprot.org/uniprotkb/P23769/entry</a> )
别名	GATA2,DCML,IMD21,MONOMAC,NFE1B,GATA2 (Phospho-Ser401) Antibody,GATA-binding protein 2,GATA2 (Phospho-Ser401)

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)