

— ABLYBIO, Help Your Research



GATA3 (Phospho-Ser308) Antibody

货号: **AYP4158**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC ELISA
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	47kDa/48kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T
细胞定位	Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human GATA3 (Phospho-Ser308).
------	----------------------------------------------------------------

靶点信息

研究背景	This gene encodes a protein which belongs to the GATA family of transcription factors. The protein contains two GATA-type zinc fingers and is an important regulator of T-cell development and plays an important role in endothelial cell biology. Defects in this gene are the cause of hypoparathyroidism with sensorineural deafness and renal dysplasia.
基因ID	2625
基因名	GATA3
Swiss	P23771 (https://www.uniprot.org/uniprotkb/P23771/entry)
别名	GATA3,HDR,HDRS,GATA3 (Phospho-Ser308) Antibody,GATA-binding factor 3,GATA3 (Phospho-Ser308)

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)