

— ABLYBIO, Help Your Research



Integrin beta 4(Phospho-Tyr1510) Antibody

货号: **AYP4190**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	106kDa/194kDa/195kDa/200kDa/202kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A431,BxPC-3,mouse pancreas
细胞定位	Cell junction,Cell membrane,Lipid-anchor,Single-pass type I membrane protein,hemidesmosome
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human Integrin beta 4(Phospho-Tyr1510).
------	--

靶点信息

研究背景	<p>Integrins are heterodimers comprised of alpha and beta subunits, that are noncovalently associated trans membrane glycoprotein receptors. Different combinations of alpha and beta polypeptides form complexes that vary in their ligand-binding specificities. Integrins mediate cell-matrix or cell-cell adhesion, and transduced signals that regulate gene expression and cell growth. This gene encodes the integrin beta 4 subunit, a receptor for the laminins. This subunit tends to associate with alpha 6 subunit and is likely to play a pivotal role in the biology of invasive carcinoma. Mutations in this gene are associated with epidermolysis bullosa with pyloric atresia. Multiple alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.</p>
基因ID	3691
基因名	ITGB4
Swiss	P16144 (https://www.uniprot.org/uniprotkb/P16144/entry)
别名	ITGB4,CD104,GP150,Integrin beta 4(Phospho-Tyr1510),Integrin beta 4(Phospho-Tyr1510) Antibody

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)