

— ABLYBIO, Help Your Research



MKP1 (Phospho-Ser359) Antibody

货号: **AYP4198**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	39kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat brain
细胞定位	Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human MKP1 (Phospho-Ser359).
------	---

靶点信息

研究背景	The expression of DUSP1 gene is induced in human skin fibroblasts by oxidative/heat stress and growth factors. It specifies a protein with structural features similar to members of the non-receptor-type protein-tyrosine phosphatase family, and which has significant amino-acid sequence similarity to a Tyr/Ser-protein phosphatase encoded by the late gene H1 of vaccinia virus. The bacterially expressed and purified DUSP1 protein has intrinsic phosphatase activity, and specifically inactivates mitogen-activated protein (MAP) kinase in vitro by the concomitant dephosphorylation of both its phosphothreonine and phosphotyrosine residues. Furthermore, it suppresses the activation of MAP kinase by oncogenic ras in extracts of <i>Xenopus</i> oocytes. Thus, DUSP1 may play an important role in the human cellular response to environmental stress as well as in the negative regulation of cellular proliferation.
基因ID	1843
基因名	DUSP1
Swiss	P28562 (https://www.uniprot.org/uniprotkb/P28562/entry)
别名	DUSP1,CL100,HVH1,MKP-1,MKP1,PTPN10,MKP1 (Phospho-Ser359) Antibody,Dual specificity protein phosphatase hVH1,Mitogen-activated protein kinase phosphatase 1,Protein-tyrosine phosphatase CL100,MKP1 (Phospho-Ser359)

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)