

— ABLYBIO, Help Your Research



# TIS11B (Phospho-Ser92) Antibody

货号: **AYP4226**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	36kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SKOV3,Rat liver
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human TIS11B (Phospho-Ser92).
------	--

## 靶点信息

研究背景	This gene is a member of the TIS11 family of early response genes, which are induced by various agonists such as the phorbol ester TPA and the polypeptide mitogen EGF. This gene is well conserved across species and has a promoter that contains motifs seen in other early-response genes. The encoded protein contains a distinguishing putative zinc finger domain with a repeating cys-his motif. This putative nuclear transcription factor most likely functions in regulating the response to growth factors. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
基因ID	677
基因名	ZFP36L1
Swiss	Q07352 ( <a href="https://www.uniprot.org/uniprotkb/Q07352/entry">https://www.uniprot.org/uniprotkb/Q07352/entry</a> )
别名	ZFP36L1, BRF1, Berg36, ERF-1, ERF1, RNF162B, TIS11B, cMG1, TIS11B (Phospho-Ser92) Antibody, Butyrate response factor 1, EGF-response factor 1, TPA-induced sequence 11b, Zinc finger protein 36, C3H1 type-like 1, TIS11B (Phospho-Ser92)

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)