

— ABLYBIO, Help Your Research



# SNAI1 (Phospho-Ser246) Antibody

货号: AYP4235

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200 <b>IF/ICC:</b> 1:50 - 1:200
理论分子量	29kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293F
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human SNAI1 (Phospho-Ser246).
------	--

## 靶点信息

研究背景	The Drosophila embryonic protein snail is a zinc finger transcriptional repressor which downregulates the expression of ectodermal genes within the mesoderm. The nuclear protein encoded by this gene is structurally similar to the Drosophila snail protein, and is also thought to be critical for mesoderm formation in the developing embryo. At least two variants of a similar processed pseudogene have been found on chromosome 2.
基因ID	6615
基因名	SNAI1
Swiss	O95863 ( <a href="https://www.uniprot.org/uniprotkb/O95863/entry">https://www.uniprot.org/uniprotkb/O95863/entry</a> )
别名	SNAI1,SLUGH2,SNA,SNAH,SNAIL,SNAIL1,dj710H13.1,SNAI1 (Phospho-Ser246) Antibody,Protein snail homolog 1,SNAI1 (Phospho-Ser246)

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)