

— ABLYBIO, Help Your Research



p62 Dok (Phospho-Tyr362) Antibody

货号: **AYP4270**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	19kDa/37kDa/52kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,U-937,Jurkat,Mouse lung,Mouse thymus
细胞定位	Cytoplasm,Nucleus,perinuclear region
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human p62 Dok (Phospho-Tyr362).
------	--

靶点信息

研究背景	The protein encoded by this gene is part of a signal transduction pathway downstream of receptor tyrosin e kinases. The encoded protein is a scaffold protein that helps form a platform for the assembly of multip rotein signaling complexes. Several transcript variants encoding different isoforms have been found for th is gene.
基因ID	1796
基因名	DOK1
Swiss	Q99704 (https://www.uniprot.org/uniprotkb/Q99704/entry)
别名	DOK1,P62DOK,pp62,p62 Dok (Phospho-Tyr362) Antibody,Downstream of tyrosine kinase 1,p62(dok),p62 Dok (Phospho-Tyr362)

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)