

— ABLYBIO, Help Your Research



FOXO1A (Phospho-Ser256) Antibody

货号: **AYP4291**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	69kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver,Mouse brain,Rat kidney
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human FOXO1A (Phospho-Ser256).
------	---

靶点信息

研究背景	This gene belongs to the forkhead family of transcription factors which are characterized by a distinct forkhead domain. The specific function of this gene has not yet been determined; however, it may play a role in myogenic growth and differentiation. Translocation of this gene with PAX3 has been associated with alveolar rhabdomyosarcoma.
基因ID	2308
基因名	FOXO1
Swiss	Q12778 (https://www.uniprot.org/uniprotkb/Q12778/entry)
别名	FKH1,FKHR,FOXO1A,FOXO1,FOXO1A (Phospho-Ser256) Antibody,Forkhead box protein O1A,Forkhead in rhabdomyosarcoma,FOXO1A (Phospho-Ser256)

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)