

— ABLYBIO, Help Your Research



FOXO4 (Phospho-Ser197) Antibody

货号: **AYP4293**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	47kDa/53kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293F,Jurkat
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human FOXO4 (Phospho-Ser197).
------	--

靶点信息

研究背景	This gene encodes a member of the O class of winged helix/forkhead transcription factor family. Proteins encoded by this class are regulated by factors involved in growth and differentiation indicating they play a role in these processes. A translocation involving this gene on chromosome X and the homolog of the Drosophila trithorax gene, encoding a DNA binding protein, located on chromosome 11 is associated with leukemia. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	4303
基因名	FOXO4
Swiss	P98177 (https://www.uniprot.org/uniprotkb/P98177/entry)
别名	AFX,AFX1,MLLT7,FOXO4,FOXO4 (Phospho-Ser197) Antibody,Fork head domain transcription factor AFX1,FOXO4 (Phospho-Ser197)

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)