

— ABLYBIO, Help Your Research



# GluR1 (Phospho-Ser849) Antibody

货号: **AYP4298**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200
理论分子量	92kDa/93kDa/101kDa/102kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Rat brain
细胞定位	Cell junction,Cell membrane,Cell projection,Endoplasmic reticulum membrane,Multi-pass membrane protein,dendrite,dendritic spine,postsynaptic cell membrane,postsynaptic density,synapse
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human GluR1 (Phospho-Ser849).
------	--

## 靶点信息

研究背景	Glutamate receptors are the predominant excitatory neurotransmitter receptors in the mammalian brain and are activated in a variety of normal neurophysiologic processes. These receptors are heteromeric protein complexes with multiple subunits, each possessing transmembrane regions, and all arranged to form a ligand-gated ion channel. The classification of glutamate receptors is based on their activation by different pharmacologic agonists. This gene belongs to a family of alpha-amino-3-hydroxy-5-methyl-4-isoxazole propionate (AMPA) receptors. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
基因ID	2890
基因名	GRIA1
Swiss	P42261 ( <a href="https://www.uniprot.org/uniprotkb/P42261/entry">https://www.uniprot.org/uniprotkb/P42261/entry</a> )
别名	GRIA1, GLUH1, GLUR1, GLURA, GluA1, HBGR1, GluR1 (Phospho-Ser849) Antibody, AMPA-selective glutamate receptor 1, GluR-A, GluR-K1, Glutamate receptor ionotropic, AMPA 1, GluR1 (Phospho-Ser849)

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)