

— ABLYBIO, Help Your Research



GSK3 α (Phospho-Ser21) Antibody

货号: **AYP4299**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IP ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	50kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,A-431,HeLa,B-cell,Mouse liver
细胞定位	apical dendrite,axon,beta-catenin destruction complex,cytoplasm,cytosol,mitochondrion,nucleus,postsynapse,proximal dendrite
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human GSK3 α (Phospho-Ser21).
------	---

靶点信息

研究背景	This gene encodes a multifunctional Ser/Thr protein kinase that is implicated in the control of several regulatory proteins including glycogen synthase, and transcription factors, such as JUN. It also plays a role in the WNT and PI3K signaling pathways, as well as regulates the production of beta-amyloid peptides associated with Alzheimer's disease.
基因ID	2931
基因名	GSK3A
Swiss	P49840 (https://www.uniprot.org/uniprotkb/P49840/entry)
别名	GSK3 alpha,GSK3A,GSK3 α (Phospho-Ser21) Antibody,Serine/threonine-protein kinase GSK3A,GSK3 α (Phospho-Ser21)

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)