

— ABLYBIO, Help Your Research



HSP27 (Phospho-Ser82) Antibody

货号: **AYP4310**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	22kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	Cytoplasm,Nucleus,cytoskeleton,spindle
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human HSP27 (Phospho-Ser82).
------	---

靶点信息

研究背景	The protein encoded by this gene is induced by environmental stress and developmental changes. The encoded protein is involved in stress resistance and actin organization and translocates from the cytoplasm to the nucleus upon stress induction. Defects in this gene are a cause of Charcot-Marie-Tooth disease type 2F (CMT2F) and distal hereditary motor neuropathy (dHMN).
基因ID	3315
基因名	HSPB1
Swiss	P04792 (https://www.uniprot.org/uniprotkb/P04792/entry)
别名	HSPB1,CMT2F,HEL-S-102,HMN2B,HS.76067,HSP27,HSP28,Hsp25,SRP27,HSP27 (Phospho-Ser82) Antibody,28 kDa heat shock protein,Estrogen-regulated 24 kDa protein,Heat shock 27 kDa protein,Heat shock protein family B member 1,Stress-responsive protein 27

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)