

— ABLYBIO, Help Your Research



IGF1R (Phospho-Tyr1165/Tyr1166) Antibody

货号: **AYP4312**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC ELISA
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	154kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,HeLa,A375
细胞定位	Cell membrane,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human IGF1R (Phospho-Tyr1165/Tyr1166).
------	---

靶点信息

研究背景	This receptor binds insulin-like growth factor with a high affinity. It has tyrosine kinase activity. The insulin-like growth factor I receptor plays a critical role in transformation events. Cleavage of the precursor gene rates alpha and beta subunits. It is highly overexpressed in most malignant tissues where it functions as an anti-apoptotic agent by enhancing cell survival. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.
基因ID	3480
基因名	IGF1R
Swiss	P08069 (https://www.uniprot.org/uniprotkb/P08069/entry)
别名	IGF1R,CD221,IGFIR,IGFR,JTK13,IGF1R (Phospho-Tyr1165/Tyr1166) Antibody,Insulin-like growth factor I receptor,IGF1R (Phospho-Tyr1165/Tyr1166)

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)