

— ABLYBIO, Help Your Research



# IRS1 (Phospho-Ser639) Antibody

货号: **AYP4318**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200
理论分子量	131kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,293T
细胞定位	caveola,cytoplasm,cytosol,nucleoplasm,nucleus,plasma membrane
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human IRS1 (Phospho-Ser639).
------	---

## 靶点信息

研究背景	This gene encodes a protein which is phosphorylated by insulin receptor tyrosine kinase. Mutations in this gene are associated with type II diabetes and susceptibility to insulin resistance.
------	--

基因ID	3667
基因名	IRS1
Swiss	P35568 ( <a href="https://www.uniprot.org/uniprotkb/P35568/entry">https://www.uniprot.org/uniprotkb/P35568/entry</a> )
别名	HIRS-1,IRS1,IRS1 (Phospho-Ser639),IRS1 (Phospho-Ser639) Antibody

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)