

— ABLYBIO, Help Your Research



# MEF2A (Phospho-Thr319) Antibody

货号: **AYP4332**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200 <b>IF/ICC:</b> 1:50 - 1:200
理论分子量	45kDa/46kDa/53kDa/54kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Raji,MCF7,Mouse spleen,C6
细胞定位	Nucleus
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human MEF2A (Phospho-Thr319).
------	--

## 靶点信息

研究背景	The protein encoded by this gene is a DNA-binding transcription factor that activates many muscle-specific, growth factor-induced, and stress-induced genes. The encoded protein can act as a homodimer or as a heterodimer and is involved in several cellular processes, including muscle development, neuronal differentiation, cell growth control, and apoptosis. Defects in this gene could be a cause of autosomal dominant coronary artery disease 1 with myocardial infarction (ADCAD1). Several transcript variants encoding different isoforms have been found for this gene.
基因ID	4205
基因名	MEF2A
Swiss	Q02078 ( <a href="https://www.uniprot.org/uniprotkb/Q02078/entry">https://www.uniprot.org/uniprotkb/Q02078/entry</a> )
别名	MEF2A,ADCAD1,RSRFC4,RSRFC9,mef2,MEF2A (Phospho-Thr319) Antibody,Serum response factor-like protein 1,MEF2A (Phospho-Thr319)

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))