

— ABLYBIO, Help Your Research



Paxillin (Phospho-Tyr31) Antibody

货号: **AYP4356**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	46kDa/60kDa/64kDa/66kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,3T3
细胞定位	Cell junction,Cytoplasm,cell cortex,cytoskeleton,focal adhesion
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human Paxillin (Phospho-Tyr31).
------	--

靶点信息

研究背景	This gene encodes a cytoskeletal protein involved in actin-membrane attachment at sites of cell adhesion to the extracellular matrix (focal adhesion). Alternatively spliced transcript variants encoding different isoforms have been described for this gene. These isoforms exhibit different expression pattern, and have different biochemical, as well as physiological properties (PMID:9054445).
基因ID	5829
基因名	PXN
Swiss	P49023 (https://www.uniprot.org/uniprotkb/P49023/entry)
别名	Paxillin,PXN,Paxillin (Phospho-Tyr31),Paxillin (Phospho-Tyr31) Antibody

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)