

— ABLYBIO, Help Your Research



PDK1 (Phospho-Ser241) Antibody

货号: **AYP4358**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	48kDa/50kDa/58kDa/60kDa/63kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HCT116,Mouse brain,Rat brain
细胞定位	Cell junction,Cell membrane,Cytoplasm,Nucleus,Peripheral membrane protein,focal adhesion
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human PDK1 (Phospho-Ser241).
------	---

靶点信息

研究背景	PDPK1 (3-phosphoinositide dependent protein kinase 1) is a protein-coding gene. GO annotations related to this gene include protein serine/threonine kinase activity and insulin receptor binding.
------	--

基因ID	5170
基因名	PDPK1
Swiss	O15530 (https://www.uniprot.org/uniprotkb/O15530/entry)
别名	PDPK1,PDK1,PDPK2,PDPK2P,PRO0461,PDK1 (Phospho-Ser241),PDK1 (Phospho-Ser241) Antibody

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)