

— ABLYBIO, Help Your Research



# PKC theta (Phospho-Ser676) Antibody

货号: **AYP4359**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200
理论分子量	67kDa/74kDa/81kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,Rat testis
细胞定位	Cell membrane,Cytoplasm,Peripheral membrane protein
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human PKC theta (Phospho-Ser676).
------	--

## 靶点信息

研究背景	Protein kinase C (PKC) is a family of serine- and threonine-specific protein kinases that can be activated by calcium and the second messenger diacylglycerol. PKC family members phosphorylate a wide variety of protein targets and are known to be involved in diverse cellular signaling pathways. PKC family members also serve as major receptors for phorbol esters, a class of tumor promoters. Each member of the PKC family has a specific expression profile and is believed to play a distinct role. The protein encoded by this gene is one of the PKC family members. It is a calcium-independent and phospholipid-dependent protein kinase. This kinase is important for T-cell activation. It is required for the activation of the transcription factors NF-kappaB and AP-1, and may link the T cell receptor (TCR) signaling complex to the activation of the transcription factors.
基因ID	5588
基因名	PRKCQ
Swiss	Q04759 ( <a href="https://www.uniprot.org/uniprotkb/Q04759/entry">https://www.uniprot.org/uniprotkb/Q04759/entry</a> )
别名	PRKCQ,PRKCT,nPKC-theta,PKC theta (Phospho-Ser676),PKC theta (Phospho-Ser676) Antibody

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)