

— ABLYBIO, Help Your Research



Src (Phospho-Tyr418) Antibody

货号: **AYP4373**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	59kDa/60kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SH-SY5Y
细胞定位	Cell membrane,Cytoplasm,Mitochondrion inner membrane,Nucleus,cytoskeleton
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human Src (Phospho-Tyr418).
------	--

靶点信息

研究背景	This gene is highly similar to the v-src gene of Rous sarcoma virus. This proto-oncogene may play a role in the regulation of embryonic development and cell growth. The protein encoded by this gene is a tyrosine-protein kinase whose activity can be inhibited by phosphorylation by c-SRC kinase. Mutations in this gene could be involved in the malignant progression of colon cancer. Two transcript variants encoding the same protein have been found for this gene.
基因ID	6714
基因名	SRC
Swiss	P12931 (https://www.uniprot.org/uniprotkb/P12931/entry)
别名	ASV,SRC1,THC6,c-SRC,p60-Src,SRC,Src (Phospho-Tyr418) Antibody,Proto-oncogene c-Src,pp60c-src,Src (Phospho-Tyr418)

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)