

— ABLYBIO, Help Your Research



# Stathmin 1 (Phospho-Ser25) Antibody

货号: **AYP4381**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200 <b>IF/ICC:</b> 1:50 - 1:200
理论分子量	17kDa/19kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,C6
细胞定位	Cytoplasm,cytoskeleton
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human Stathmin 1 (Phospho-Ser25).
------	--

## 靶点信息

研究背景	This gene belongs to the stathmin family of genes. It encodes a ubiquitous cytosolic phosphoprotein proposed to function as an intracellular relay integrating regulatory signals of the cellular environment. The encoded protein is involved in the regulation of the microtubule filament system by destabilizing microtubules. It prevents assembly and promotes disassembly of microtubules. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	3925
基因名	STMN1
Swiss	P16949 ( <a href="https://www.uniprot.org/uniprotkb/P16949/entry">https://www.uniprot.org/uniprotkb/P16949/entry</a> )
别名	STMN1,C1orf215,LAP18,Lag,OP18,PP17,PP19,PR22,SMN,stathmin,Stathmin 1 (Phospho-Ser25) Antibody,L eukemia-associated phosphoprotein p18,Metablastin,Oncoprotein 18,Phosphoprotein p19,Prosolin,Protein Pr22,Stathmin 1 (Phospho-Ser25)

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)