

Stathmin 1 (Phospho-Ser38) Antibody

货号: **AYP4382**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	17kDa/19kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,C6
细胞定位	Cytoplasm,cytoskeleton
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human Stathmin 1 (Phospho-Ser38).
------	--

靶点信息

研究背景	This gene belongs to the stathmin family of genes. It encodes a ubiquitous cytosolic phosphoprotein proposed to function as an intracellular relay integrating regulatory signals of the cellular environment. The encoded protein is involved in the regulation of the microtubule filament system by destabilizing microtubules. It prevents assembly and promotes disassembly of microtubules. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	3925
基因名	STMN1

Swiss	P16949
别名	STMN1;C1orf215;LAP18;Lag;OP18;PP17;PP19;PR22;SMN;stathmin

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn