

— ABLYBIO, Help Your Research



# VEGFR2 (Phospho-Tyr1175) Antibody

货号: **AYP4399**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200
理论分子量	75kDa/79kDa/151kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart
细胞定位	Cell junction,Cell membrane,Cytoplasm,Cytoplasmic vesicle,Early endosome,Endoplasmic reticulum,Nucleus,Secreted,Single-pass type I membrane protein
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human VEGFR2 (Phospho-Tyr1175).
------	--

## 靶点信息

研究背景	Vascular endothelial growth factor (VEGF) is a major growth factor for endothelial cells. This gene encodes one of the two receptors of the VEGF. This receptor, known as kinase insert domain receptor, is a type II I receptor tyrosine kinase. It functions as the main mediator of VEGF-induced endothelial proliferation, survival, migration, tubular morphogenesis and sprouting. The signalling and trafficking of this receptor are regulated by multiple factors, including Rab GTPase, P2Y purine nucleotide receptor, integrin alphaVbeta3, T-cell protein tyrosine phosphatase, etc.. Mutations of this gene are implicated in infantile capillary hemangiomas.
基因ID	3791
基因名	KDR
Swiss	P35968 ( <a href="https://www.uniprot.org/uniprotkb/P35968/entry">https://www.uniprot.org/uniprotkb/P35968/entry</a> )
别名	CD309,FLK1,VEGFR,VEGFR2,VEGF Receptor 2,KDR,VEGFR2 (Phospho-Tyr1175) Antibody,Fetal liver kinase 1,Kinase insert domain receptor,Protein-tyrosine kinase receptor flk-1,VEGFR2 (Phospho-Tyr1175)

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)