

— ABLYBIO, Help Your Research



# ZAP70 (Phospho-Tyr319) Antibody

货号: **AYP4401**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200 <b>IF/ICC:</b> 1:50 - 1:200
理论分子量	35kDa/55kDa/69kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,Mouse thymus
细胞定位	Cell membrane,Cytoplasm,Peripheral membrane protein
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human ZAP70 (Phospho-Tyr319).
------	--

## 靶点信息

研究背景	This gene encodes an enzyme belonging to the protein tyrosine kinase family, and it plays a role in T-cell development and lymphocyte activation. This enzyme, which is phosphorylated on tyrosine residues upon T-cell antigen receptor (TCR) stimulation, functions in the initial step of TCR-mediated signal transduction in combination with the Src family kinases, Lck and Fyn. This enzyme is also essential for thymocyte development. Mutations in this gene cause selective T-cell defect, a severe combined immunodeficiency disease characterized by a selective absence of CD8-positive T-cells. Two transcript variants that encode different isoforms have been found for this gene.
基因ID	7535
基因名	ZAP70
Swiss	P43403 ( <a href="https://www.uniprot.org/uniprotkb/P43403/entry">https://www.uniprot.org/uniprotkb/P43403/entry</a> )
别名	ZAP70,ADMIO2,IMD48,SRK,STCD,STD,TZK,ZAP-70,ZAP70 (Phospho-Tyr319) Antibody,70 kDa zeta-chain associated protein,Syk-related tyrosine kinase,ZAP70 (Phospho-Tyr319)

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)