

— ABLYBIO, Help Your Research



# Ezrin (Phospho-Tyr146) Antibody

货号: **AYP4408**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	69kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MDA
细胞定位	Apical cell membrane,Cell projection,Cytoplasm,Cytoplasmic side,Peripheral membrane protein,cell cortex,cytoskeleton,microvillus membrane,ruffle membrane
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human Ezrin (Phospho-Tyr146).
------	--

## 靶点信息

研究背景	The cytoplasmic peripheral membrane protein encoded by this gene functions as a protein-tyrosine kinase substrate in microvilli. As a member of the ERM protein family, this protein serves as an intermediate between the plasma membrane and the actin cytoskeleton. This protein plays a key role in cell surface structure adhesion, migration and organization, and it has been implicated in various human cancers. A pseudogene located on chromosome 3 has been identified for this gene. Alternatively spliced variants have also been described for this gene.
基因ID	7430
基因名	EZR
Swiss	P15311 ( <a href="https://www.uniprot.org/uniprotkb/P15311/entry">https://www.uniprot.org/uniprotkb/P15311/entry</a> )
别名	CVIL,CVL,HEL-S-105,VIL2,EZR,Ezrin,Ezrin (Phospho-Tyr146) Antibody,Cytovillin,Villin-2,p81,Ezrin (Phospho-Tyr146)

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))