

— ABLYBIO, Help Your Research



IKK- γ (Phospho-Ser376) Antibody

货号: **AYP4446**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	36kDa/48kDa/55kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,Mouse liver
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human IKK- γ (Phospho-Ser376).
------	--

靶点信息

研究背景	This gene encodes the regulatory subunit of the inhibitor of kappaB kinase (IKK) complex, which activates NF-kappaB resulting in activation of genes involved in inflammation, immunity, cell survival, and other pathways. Mutations in this gene result in incontinentia pigmenti, hypohidrotic ectodermal dysplasia, and several other types of immunodeficiencies. A pseudogene highly similar to this locus is located in an adjacent region of the X chromosome.
基因ID	8517
基因名	IKBKG
Swiss	Q9Y6K9 (https://www.uniprot.org/uniprotkb/Q9Y6K9/entry)
别名	IKBKG,AMCBX1,FIP-3,FIP3,Fip3p,IKK-gamma,IKKAP1,IKKG,IMD33,IP,IP1,IP2,IPD2,NEMO,ZC2HC9,IKK- γ (Phospho-Ser376) Antibody,IkB kinase-associated protein 1,Inhibitor of nuclear factor kappa-B kinase subunit gamma,NF-kappa-B essential modifier

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)