

— ABLYBIO, Help Your Research



# IRAK1 (Phospho-Ser376) Antibody

货号: **AYP4451**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200
理论分子量	68kDa/73kDa/74kDa/76kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,A-549,Mouse liver,Mouse kidney,Mouse skeletal muscle,Rat liver,Rat kidney,Rat skeletal muscle
细胞定位	Cytoplasm,Lipid droplet,Nucleus
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human IRAK1 (Phospho-Ser376).
------	--

## 靶点信息

研究背景	This gene encodes the interleukin-1 receptor-associated kinase 1, one of two putative serine/threonine kinases that become associated with the interleukin-1 receptor (IL1R) upon stimulation. This gene is partially responsible for IL1-induced upregulation of the transcription factor NF-kappa B. Alternatively spliced transcripts encoding different isoforms have been found for this gene.
基因ID	3654
基因名	IRAK1
Swiss	P51617 ( <a href="https://www.uniprot.org/uniprotkb/P51617/entry">https://www.uniprot.org/uniprotkb/P51617/entry</a> )
别名	IRAK1,IRAK,pelle,IRAK1 (Phospho-Ser376),IRAK1 (Phospho-Ser376) Antibody

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)