

# JAK2 (Phospho-Tyr931) Antibody

货号: **AYP4452**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200 <b>IF/ICC:</b> 1:50 - 1:200
理论分子量	130kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT29
细胞定位	Cytoplasm,Endomembrane system,Nucleus,Peripheral membrane protein
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human JAK2 (Phospho-Tyr931).
------	---

## 靶点信息

研究背景	This gene product is a protein tyrosine kinase involved in a specific subset of cytokine receptor signaling pathways. It has been found to be constitutively associated with the prolactin receptor and is required for responses to gamma interferon. Mice that do not express an active protein for this gene exhibit embryonic lethality associated with the absence of definitive erythropoiesis.
基因ID	3717
基因名	JAK2

Swiss	O60674
别名	JAK2;JTK10;THCYT3

产品验证

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)