

— ABLYBIO, Help Your Research



STK39 (Phospho-Ser325) Antibody

货号: **AYP4482**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	68kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,SH-SY5Y,Rat testis,Mouse testis
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human STK39 (Phospho-Ser325).
------	--

靶点信息

研究背景	This gene encodes a serine/threonine kinase that is thought to function in the cellular stress response pathway. The kinase is activated in response to hypotonic stress, leading to phosphorylation of several cation-chloride-coupled cotransporters. The catalytically active kinase specifically activates the p38 MAP kinase pathway, and its interaction with p38 decreases upon cellular stress, suggesting that this kinase may serve as an intermediate in the response to cellular stress. [provided by RefSeq, Jul 2008]
基因ID	27347
基因名	STK39
Swiss	Q9UEW8 (https://www.uniprot.org/uniprotkb/Q9UEW8/entry)
别名	DCHT,PASK,SPAK,STK39 (Phospho-Ser325) Antibody,STK39,Serine/threonine-protein kinase 39,STK39 (Phospho-Ser325)

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)