

— ABLYBIO, Help Your Research



Trk B (Phospho-Tyr706+Tyr707) Antibody

货号: **AYP4486**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	35kDa/53kDa/59kDa/60kDa/81kDa/91kDa/93kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Rat brain
细胞定位	Cell membrane,Endosome membrane,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human Trk B (Phospho-Tyr706+Tyr707).
------	---

靶点信息

研究背景	This gene encodes a member of the neurotrophic tyrosine receptor kinase (NTRK) family. This kinase is a membrane-bound receptor that, upon neurotrophin binding, phosphorylates itself and members of the M APK pathway. Signalling through this kinase leads to cell differentiation. Mutations in this gene have been associated with obesity and mood disorders. Alternative splicing results in multiple transcript variants.
基因ID	4915
基因名	NTRK2
Swiss	Q16620 (https://www.uniprot.org/uniprotkb/Q16620/entry)
别名	NTRK2,GP145-TrkB,TRKB,tk-B,OBHD,Trk B (Phospho-Tyr706+Tyr707) Antibody,Neurotrophic tyrosine kinase receptor type 2,TrkB tyrosine kinase,Tropomyosin-related kinase B,Trk B (Phospho-Tyr706+Tyr707)

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)