

— ABLYBIO, Help Your Research



iNOS (Phospho-Tyr151) Antibody

货号: **AYP4502**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	131kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	NIH/3T3,C6,C2C12
细胞定位	cortical cytoskeleton,cytoplasm,cytosol,nucleoplasm,nucleus,perinuclear region of cytoplasm,peroxisoma l matrix,peroxisome,plasma membrane
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human iNOS (Phospho-Tyr151).
------	---

靶点信息

研究背景	Nitric oxide is a reactive free radical which acts as a biologic mediator in several processes, including neurotransmission and antimicrobial and antitumoral activities. This gene encodes a nitric oxide synthase which is expressed in liver and is inducible by a combination of lipopolysaccharide and certain cytokines. Three related pseudogenes are located within the Smith-Magenis syndrome region on chromosome 17.
基因ID	4843
基因名	NOS2
Swiss	P35228 (https://www.uniprot.org/uniprotkb/P35228/entry)
别名	NOS, INOS, NOS2A, HEP-NOS, NOS2, iNOS (Phospho-Tyr151) Antibody, Hepatocyte NOS, Inducible NO synthase, NOS type II, Peptidyl-cysteine S-nitrosylase NOS2, iNOS (Phospho-Tyr151)

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)