

— ABLYBIO, Help Your Research



Connexin 43 (Phospho-Tyr265) Antibody

货号: **AYP4503**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	43kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Mouse heart,Rat heart
细胞定位	Cell junction,Cell membrane,Endoplasmic reticulum,Multi-pass membrane protein,gap junction
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human Connexin 43 (Phospho-Tyr265).
------	--

靶点信息

研究背景	This gene is a member of the connexin gene family. The encoded protein is a component of gap junctions, which are composed of arrays of intercellular channels that provide a route for the diffusion of low molecular weight materials from cell to cell. The encoded protein is the major protein of gap junctions in the heart that are thought to have a crucial role in the synchronized contraction of the heart and in embryonic development. A related intronless pseudogene has been mapped to chromosome 5. Mutations in this gene have been associated with oculodentodigital dysplasia, autosomal recessive craniometaphyseal dysplasia and heart malformations.
基因ID	2697
基因名	GJA1
Swiss	P17302 (https://www.uniprot.org/uniprotkb/P17302/entry)
别名	GJA1,AVSD3,CMDR,CX43,EKVP,GJAL,HLHS1,HSS,ODDD,PPKCA,EKVP3,Connexin 43 (Phospho-Tyr265) Antibody,Connexin-43,Gap junction 43 kDa heart protein,Connexin 43 (Phospho-Tyr265)

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)