

— ABLYBIO, Help Your Research



ITCH (Phospho-Tyr420) Antibody

货号: **AYP4513**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	103kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,THP-1,NIH/3T3,C6,Mouse lung,Mouse kidney,Rat lung,Rat spleen
细胞定位	cell cortex,cytoplasm,cytoplasmic vesicle,cytosol,extracellular exosome,nucleoplasm,plasma membrane
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human ITCH (Phospho-Tyr420).
------	---

靶点信息

研究背景	This gene encodes a member of the Nedd4 family of HECT domain E3 ubiquitin ligases. HECT domain E3 ubiquitin ligases transfer ubiquitin from E2 ubiquitin-conjugating enzymes to protein substrates, thus targeting specific proteins for lysosomal degradation. The encoded protein plays a role in multiple cellular processes including erythroid and lymphoid cell differentiation and the regulation of immune responses. Mutations in this gene are a cause of syndromic multisystem autoimmune disease. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
基因ID	83737
基因名	ITCH
Swiss	Q96J02 (https://www.uniprot.org/uniprotkb/Q96J02/entry)
别名	ITCH,ADMFD,AIF4,AIP4,NAPP1,ITCH (Phospho-Tyr420) Antibody,Atrophin-1-interacting protein 4,HECT-type E3 ubiquitin transferase Itchy homolog,NFE2-associated polypeptide 1,ITCH (Phospho-Tyr420)

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)