

— ABLYBIO, Help Your Research



Caveolin 2 (Phospho-Tyr27) Antibody

货号: **AYP4517**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	18kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Mouse lung
细胞定位	Cell membrane,Cytoplasm,Golgi apparatus membrane,Membrane,Nucleus,Peripheral membrane protein, caveola
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human Caveolin 2 (Phospho-Tyr27).
------	--

靶点信息

研究背景	The protein encoded by this gene is a major component of the inner surface of caveolae, small invaginations of the plasma membrane, and is involved in essential cellular functions, including signal transduction, lipid metabolism, cellular growth control and apoptosis. This protein may function as a tumor suppressor. This gene and related family member (CAV1) are located next to each other on chromosome 7, and express colocalizing proteins that form a stable hetero-oligomeric complex. Alternatively spliced transcript variants encoding different isoforms have been identified for this gene. Additional isoforms resulting from the use of alternate in-frame translation initiation codons have also been described, and shown to have preferential localization in the cell (PMID:11238462). [provided by RefSeq, May 2011]
基因ID	858
基因名	CAV2
Swiss	P51636 (https://www.uniprot.org/uniprotkb/P51636/entry)
别名	CAV,Caveolin 2 (Phospho-Tyr27),Caveolin 2 (Phospho-Tyr27) Antibody,CAV2

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)