

— ABLYBIO, Help Your Research



Neuro D (Phospho-Ser274) Antibody

货号: **AYP4562**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	39kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	C6,SH-SY5Y
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human Neuro D (Phospho-Ser274).
------	--

靶点信息

研究背景	This gene encodes a member of the NeuroD family of basic helix-loop-helix (bHLH) transcription factors. The protein forms heterodimers with other bHLH proteins and activates transcription of genes that contain a specific DNA sequence known as the E-box. It regulates expression of the insulin gene, and mutations in this gene result in type II diabetes mellitus.
基因ID	4760
基因名	NEUROD1
Swiss	Q13562 (https://www.uniprot.org/uniprotkb/Q13562/entry)
别名	NEUROD1,BETA2,BHF-1,MODY6,NEUROD,bHLHa3,Neuro D (Phospho-Ser274) Antibody,Class A basic helix-loop-helix protein 3,Neuro D (Phospho-Ser274)

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)