

— ABLYBIO, Help Your Research



# NudC (Phospho-Ser326) Antibody

货号: **AYP4565**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200
理论分子量	42kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,293T,HepG2,A-549,Mouse lung,Mouse brain,Rat brain
细胞定位	Cytoplasm,Nucleus,cytoskeleton
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human NudC (Phospho-Ser326).
------	---------------------------------------------------------------

## 靶点信息

研究背景	This gene encodes a nuclear distribution protein that plays an essential role in mitosis and cytokinesis. The encoded protein is involved in spindle formation during mitosis and in microtubule organization during cytokinesis. Pseudogenes of this gene are found on chromosome 2. [provided by RefSeq, Feb 2012]
基因ID	10726
基因名	NUDC
Swiss	Q9Y266 ( <a href="https://www.uniprot.org/uniprotkb/Q9Y266/entry">https://www.uniprot.org/uniprotkb/Q9Y266/entry</a> )
别名	HNUDC,MNUDC,NPD011,NudC (Phospho-Ser326) Antibody,NUDC,Nuclear distribution protein C homolog, NudC (Phospho-Ser326)

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)