

— ABLYBIO, Help Your Research



# p47 phox (Phospho-Ser345) Antibody

货号: **AYP4567**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200
理论分子量	30kDa/44kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Raji,Daudi,Jurkat(Negative control)
细胞定位	Cytoplasm,Cytoplasmic side,Membrane,Peripheral membrane protein,cytosol
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human p47 phox (Phospho-Ser345).
------	---

## 靶点信息

研究背景	The protein encoded by this gene is a 47 kDa cytosolic subunit of neutrophil NADPH oxidase. This oxidase is a multicomponent enzyme that is activated to produce superoxide anion. Mutations in this gene have been associated with chronic granulomatous disease.
基因ID	653361
基因名	NCF1
Swiss	P14598 ( <a href="https://www.uniprot.org/uniprotkb/P14598/entry">https://www.uniprot.org/uniprotkb/P14598/entry</a> )
别名	NCF1,NCF1A,NOXO2,SH3PXD1A,p47phox,p47 phox (Phospho-Ser345) Antibody,47 kDa autosomal chronic granulomatous disease protein,47 kDa neutrophil oxidase factor,NCF-47K,Neutrophil NADPH oxidase factor 1,Nox organizer 2,Nox-organizing protein 2,SH3 and PX domain-containing protein 1A

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)