

— ABLYBIO, Help Your Research



RPS6 (Phospho-Ser235+Ser236) Antibody

货号: **AYP4573**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	28kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SKOV3,Raji,HT-29,A-431,Mouse heart
细胞定位	cytoplasm,cytoplasmic ribonucleoprotein granule,cytosol,cytosolic ribosome,endoplasmic reticulum,nucleolus,nucleoplasm,nucleus,perinuclear region of cytoplasm
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human RPS6 (Phospho-Ser235+Ser236).
------	--

靶点信息

研究背景	Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a cytoplasmic ribosomal protein that is a component of the 40S subunit. The protein belongs to the S6E family of ribosomal proteins. It is the major substrate of protein kinases in the ribosome, with subsets of five C-terminal serine residues phosphorylated by different protein kinases. Phosphorylation is induced by a wide range of stimuli, including growth factors, tumor-promoting agents, and mitogens. Dephosphorylation occurs at growth arrest. The protein may contribute to the control of cell growth and proliferation through the selective translation of particular classes of mRNA. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed throughout the genome.
基因ID	6194
基因名	RPS6
Swiss	P62753 (https://www.uniprot.org/uniprotkb/P62753/entry)
别名	RPS6,S6,RPS6 (Phospho-Ser235+Ser236) Antibody,40S ribosomal protein S6,Phosphoprotein NP33,RPS6 (Phospho-Ser235+Ser236)

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)