

— ABLYBIO, Help Your Research



# VCP (Phospho-Ser352) Antibody

货号: **AYP4583**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	89kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,K-562,SH-SY5Y,Mouse brain,Mouse spleen,Mouse kidney,Rat liver,Rat brain
细胞定位	Cytoplasm,Endoplasmic reticulum,Nucleus,cytosol
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human VCP (Phospho-Ser352).
------	--

## 靶点信息

研究背景	The protein encoded by this gene is a member of a family that includes putative ATP-binding proteins involved in vesicle transport and fusion, 26S proteasome function, and assembly of peroxisomes. This protein, as a structural protein, is associated with clathrin, and heat-shock protein Hsc70, to form a complex. It has been implicated in a number of cellular events that are regulated during mitosis, including homotypic membrane fusion, spindle pole body function, and ubiquitin-dependent protein degradation.
基因ID	7415
基因名	VCP
Swiss	P55072
别名	VCP,ALS14,CDC48,CMT2Y,HEL-220,HEL-S-70,IBMPFD,IBMPFD1,TERA,p97,VCP (Phospho-Ser352) Antibody, 15S Mg(2+)-ATPase p97 subunit,Valosin-containing protein,VCP (Phospho-Ser352)

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)