

— ABLYBIO, Help Your Research



SYK (Phospho-Tyr525) Antibody

货号: **AYP4604**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	69kDa/72kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,Mouse spleen,Rat liver
细胞定位	Cell membrane,Cytoplasm,cytosol
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human SYK (Phospho-Tyr525).
------	--

靶点信息

研究背景	This gene encodes a member of the family of non-receptor type Tyr protein kinases. This protein is widely expressed in hematopoietic cells and is involved in coupling activated immunoreceptors to downstream signaling events that mediate diverse cellular responses, including proliferation, differentiation, and phagocytosis. It is thought to be a modulator of epithelial cell growth and a potential tumour suppressor in human breast carcinomas. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
基因ID	6850
基因名	SYK
Swiss	P43405
别名	SYK,p72-Syk,SYK (Phospho-Tyr525) Antibody,Spleen tyrosine kinase,SYK (Phospho-Tyr525)

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)