

— ABLYBIO, Help Your Research



# HSP90B Antibody

货号: **AYP4623**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC IP ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200 <b>IF/ICC:</b> 1:50 - 1:200 <b>IP:</b> 1:20 - 1:50
理论分子量	83kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	K-562,DU145,SW620,COS-7,NIH/3T3,HeLa
细胞定位	Cytoplasm,Melanosome
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human HSP90B.
------	--

## 靶点信息

研究背景	This gene encodes a member of the heat shock protein 90 family; these proteins are involved in signal transduction, protein folding and degradation and morphological evolution. This gene encodes the constitutive form of the cytosolic 90 kDa heat-shock protein and is thought to play a role in gastric apoptosis and inflammation. Alternative splicing results in multiple transcript variants. Pseudogenes have been identified on multiple chromosomes.
基因ID	3326
基因名	HSP90AB1
Swiss	P08238 ( <a href="https://www.uniprot.org/uniprotkb/P08238/entry">https://www.uniprot.org/uniprotkb/P08238/entry</a> )
别名	D6S182,HSP84,HSP90B,HSPC2,HSPCB,Hsp90 beta,HSP90AB1,HSP90B Antibody,Heat shock 84 kDa,Heat shock protein family C member 3,HSPC3

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)